



[SIOT WEBSITE](#)

[Contact us](#)

[Customer section](#)

## ONE OF THE SMALLEST THERMOSTATS YOU CAN FIND

### LABORATORY THERMOSTAT



- temperature stability +/- 0.5 °C
- net volume: 45 l
- Maximum power consumption: 130W
- Temperature range from +8°C to +40°C
- External dimensions (width/depth/height): 470/440/540 mm
- Internal dimensions ( width/depth/height): 400/350/395 mm

### LABORATORY REFRIGERATOR



- temperature stability +/- 0,5 °C
- net volume: 45 l
- temperature range 0°C to +15°C
- maximum power consumption: 130W
- External dimensions (width/depth/height): 470/440/540 mm
- Internal dimensions ( width/depth/height): 400/350/395 mm



**Petr Kocurek**

Technical support

SIOT Trade s.r.o.

Tel.: +420 736 231 783

E-mail: [p.kocurek@siot.cz](mailto:p.kocurek@siot.cz)

